

## Important Stifle Fitting Information

Please review the videos from the attached link:

<https://vimeo.com/album/1604074/video/24051296>

1. Place brace on the hock at a 90 degree angle. This needs to be a snug fit, but not too tight. Fasten the bottom strap, **not too tight**.



2. Push brace onto knee. **Push** from the knee hinges. Then **squeeze** the brace inward (medial and lateral). While squeezing, fasten the back middle velcro strap **snugly**. Only one finger should fit between the strap and skin



3. Push the thigh cuff onto the leg. Squeeze cuff inward (medial and lateral) and while squeezing fasten the top back strap **snugly**. One finger should fit between strap and skin.



4. Fasten front strap. This does not need to be pulled too tight or it will just bring the brace forward off the knee.



5. Trim velcro straps  
You can move and/or trim the pads if needed.

